



sartorius stedim  
biotech

## SARTOFLOW® Slice 200 Benchtop System Welcome to the world of downscale crossflow filtration



Get 4  
cassettes  
for free

# To optimize your membrane trials...

This package includes a SARTOFLOW® Slice 200 Benchtop System, an extra pack of three pressure transmitters and a plastic fittings kit. Additionally you receive 4 complimentary Sartoflow® Slice 200 Cassettes.

The SARTOFLOW® Slice 200 Benchtop System allows you to enter the world of crossflow filtration. This system is automated and perfectly suited for R&D, process development, pre-clinical and small pilot lots. It is easy to handle and ready to use in minutes.

The features of the SARTOFLOW® Slice 200 Benchtop System allows users to automate, optimize and document the crossflow processes:

- Reliable operation through implementation of sensors & user-definable alarms
- On-board application programs that automate and optimize filtration

- Finger-tip control via the front panel with menu-driven software
- Interface directly to PC
- Real-time verification and documentation of process parameters

## Technical Specifications

- Filter Area 200 to 400 cm<sup>2</sup>
- Vessel Volume 500 ml
- Min. Operating Volume < 20 ml
- Max. System Pressure 58 psi | 4 bar
- Luer Lock connections
- Dimensions (W + D) 60 cm + 40 cm
- RS232 Communication



...get



## 4 Sartococon® Slice 200 crossflow cassettes for free

### Membranes you always wanted

Your demand is to efficiently process highly valuable products while securing the highest yield. High-performance membranes fit into your process through:

- Separation
- Concentration
- Diafiltration
- Final Purification

The Crossflow Cassettes are designed based on application and process demands. This leads to process-oriented cassette geometry, the correct membrane material and resolute process control. This helps to minimize residue build-up and prevent blockages.

# Welcome Package | Order form

**Special offer!**  
**Only available**  
**until 30.9.10**



## SARTOFLOW® Slice 200 Benchtop System

- 17525SYS-BT1 SARTOFLOW® Slice 200 Benchtop System (120V)
- 17525SYS-BT2 SARTOFLOW® Slice 200 Benchtop System (220V)

## Spare parts

- 17525SP-10 Pack of three Pressure Transmitters
- 17525SP-02 Plastic fittings kit
- 17521---037 Reusable Pressure Sensor Kit **NEW**

## Mircofiltration

Order code	Material	Membrane area	Pore size
<input type="checkbox"/> 3081860702W--SW	Hydrosart®	200 cm <sup>2</sup>	0.20 µm
<input type="checkbox"/> 3081860602W--SW	Hydrosart®	200 cm <sup>2</sup>	0.45 µm
<input type="checkbox"/> 3081545802W--SW	PESU	200 cm <sup>2</sup>	0.10 µm

## Ultrafiltration

Order code	Material	Membrane area	Cut-off
<input type="checkbox"/> 3081441902E--SW	Hydrosart®	200 cm <sup>2</sup>	2 kDa
<input type="checkbox"/> 3081442902E--SW	Hydrosart®	200 cm <sup>2</sup>	5 kDa
<input type="checkbox"/> 3081443902E--SW	Hydrosart®	200 cm <sup>2</sup>	10 kDa
<input type="checkbox"/> 3081445902E--SW	Hydrosart®	200 cm <sup>2</sup>	30 kDa
<input type="checkbox"/> 3081446802E--SW	Hydrosart®	200 cm <sup>2</sup>	100 kDa
<input type="checkbox"/> 3081460902E--SW	PESU	200 cm <sup>2</sup>	1 kDa
<input type="checkbox"/> 3081462902E--SW	PESU	200 cm <sup>2</sup>	5 kDa
<input type="checkbox"/> 3081463402E--SW	PESU	200 cm <sup>2</sup>	8 kDa
<input type="checkbox"/> 3081463902E--SW	PESU	200 cm <sup>2</sup>	10 kDa
<input type="checkbox"/> 3081465902E--SW	PESU	200 cm <sup>2</sup>	30 kDa
<input type="checkbox"/> 3081465002E--SW	PESU	200 cm <sup>2</sup>	50 kDa
<input type="checkbox"/> 3081466802E--SW	PESU	200 cm <sup>2</sup>	100 kDa
<input type="checkbox"/> 3081467902E--SW	PESU	200 cm <sup>2</sup>	300 kDa

## I want to order further Sartocan® Slice 200 Cassettes

Order Code: \_\_\_\_\_

Please fill in amount: \_\_\_\_\_

For further information please contact your local sales representative

# Customer Details

**Please complete this form and return to**

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen, Germany

Via Fax: +49.551.308.3838

Company Name: \_\_\_\_\_

Contact Person: \_\_\_\_\_

Email address: \_\_\_\_\_

Phone number: \_\_\_\_\_

Fax number: \_\_\_\_\_

Delivery address: \_\_\_\_\_

City: \_\_\_\_\_

ZIP Code: \_\_\_\_\_

Country: \_\_\_\_\_

Date: \_\_\_\_\_

Signature: \_\_\_\_\_

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen, Germany

Phone +49.551.308.0  
Fax +49.551.308.3289

[www.sartorius-stedim.com](http://www.sartorius-stedim.com)

USA +1.800.368.7178  
UK +44.1372.737159  
France +33.442.845600  
Italy +39.055.63.40.41  
Spain +34.90.2110935  
Japan +81.3.3740.5407